

Challenge

Emerging technologies for molecular microbial ecology are providing holistic views of microbial communities. However, challenges in metabolic patterns and community diversity make it difficult to effectively capture “snapshots” of microbial communities at such studies. Given that on-location preservation of samples are difficult to achieve, nucleic acids for later extraction is often the preferred method. RNA and biostaticness are both critical for sample integrity during transport and storage. Temperature conditions during transport and storage on microbial communities [1]. Reagents used for RNA extraction have been known to be inadequate for soil samples, and low quality RNA after extraction, thus affecting downstream methods.

Solution (Literally)

As a novel reagent, the MO BIO LifeGuard™ Soil Stabilizer is designed to stabilize total microbial RNA in soils at various temperatures. LifeGuard™ is a patented formulation of nucleic acids through enzyme inhibition that provides the benefit of keeping samples biostatic [2].

Materials and Methods

Soils were used at a ratio (2.5X) of 5 ml of solution to 12.5 ml of soil (quality assessment) and various ratios for quality assessment. Soils from the Dry Valleys, Antarctica (for examining microbial communities) were stored at indicated temperatures. RNA was extracted using the MO BIO PowerSoil Total RNA Isolation Kit. Samples were analyzed using agarose gel electrophoresis. Select samples were also analyzed using RT-qPCR to evaluate extraction efficiency.

RNA Preservation Solutions

RNA extracted from soil samples using LifeGuard™, RNALater, and GITC/Solution D (5M guanidinium thiocyanate) storage at 4°C and -20°C over a four-week period. High quality RNA at all time-points and temperatures was observed.


C. K. Lee¹, S. J. Kennedy², A. Rueckert¹, V. Moroney², M. N. Brolaski²,

¹Thermophile Research Unit, University of Waikato, Hamilton, New Zealand; ²MO BIO Laboratories, 3College of Earth, Ocean, and Environment, University of Delaware, Lewes,

Figure 1

Figure 1A shows an agarose gel analysis of RNA extracted from soil samples preserved with LifeGuard™, 70% EtOH, and RNALater for two weeks at -20°C. 70% EtOH, while retaining some RNA, showed significant degradation; whereas nearly no RNA can be recovered from samples preserved with LifeGuard™ or RNALater.

Figure 1B shows spectrophotometric (i.e. NanoDrop) profiles of RNA preparations from Figure 1A, and it is apparent that those of samples preserved in RNALater were highly abnormal. The degradation of RNA in 70% EtOH-preserved samples is not reflected in their NanoDrop profiles.




Quantification of Preserved RNA

Using armored RNA [3] as the internal standard for extraction efficiency, we found that the amounts of recoverable RNA declined with time and differed between various preservation solutions and temperatures. Both RNALater and GITC failed to yield any RNA beyond two weeks.

Figure 2

Figure 2A shows the relationship between total RNA and armored RNA recovery ratios. A linear correlation was found, suggesting the amount of RNA extracted is indicative of recoverable RNA in the samples. Figure 2B shows the amounts of RNA recovered over the four-week period from various preservation solutions and temperatures.



Bacteriostatic Properties

To assess the conservation of microbial community structure in LifeGuard™-preserved samples, reverse transcription terminal restriction fragment polymorphism (RT-tRFLP) analysis using 16S rRNA gene-specific primers was performed on RNA samples extracted from preserved Dry Valley soils. The data indicated that after two weeks at -20°C, the RT-tRFLP profiles were statistically indistinguishable from those of the zero-time control samples. Storage at higher temperatures (i.e., 4°C) is possible for short periods, but LifeGuard™:Soil ratios lower than 2X should be avoided. Storage at warm temperatures (e.g., 37°C) potentially leads to inconsistency in observed microbial community structure.

Conclusion

Combining analyses of RNA quality and RT-tRFLP, findings from our study suggest nucleic acid stabilizer and can effectively preserve under appropriate storage conditions. LifeGuard™ offers a high quality total microbial RNA recovery, and is directly compatible with the LifeGuard™ Isolation Kit. Additionally, because LifeGuard™ is bactericidal, it may also have applications in

References

1. R.A. Hurt, X. Qiu, et al. Appl Environ Microbiol (2010), 76:1000-1008
2. J.P. Audia, M.C. Patton, and H.H. Winkler. Appl Environ Microbiol (2004), 70:1000-1008
3. A. Rueckert and S.C. Cary. Limnol. Oceanogr. (2004), 49:1000-1008